

Product Datasheet

vegfaa antibody (orb53331)

Description	vegfaa antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Vascular endothelial growth factor A-A protein (24-188AA)
Target	vegfaa
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O73682
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.