

## Product Datasheet

### qPCR SybrMaster highROX (orb533384)

<b>Description</b>	qPCR SybrMaster highROX is designed for quantitative real-time analysis of DNA samples. The mix cont...
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	2x conc
<b>Storage</b>	store at -20 °C
<b>Note</b>	For research use only
<b>Application notes</b>	<b>Spectroscopic Propertie:</b> $\lambda_{exc}$ 494 nm (bound to DNA), $\lambda_{em}$ 521 nm, (bound to DNA).
<b>Expiration Date</b>	6 months from date of receipt.