

Product Datasheet

Yeast DNA Preparation - Solution Kit (orb533443)

Description	Yeast DNA Preparation Kit is designed for convenient and fast isolation of genomic DNA from yeast cells. The solution based system minimizes DNA fragmentation that may be problematic in other spin-column/filtration based method. Because phenol or chloroform is not used it is safe and does not produce any harmful waste. Solution based genomic DNA purification kits guarantee minimal DNA fragmentation and yield DNA sized up to 150 kb.
Storage	store at ambient temperature
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)