

Product Datasheet

SAP-Exo Kit (orb533451)

Description	The Kit contains two hydrolytic enzymes, recombinant Shrimp Alka-line Phosphatase (rSAP) and Exonucl...
Storage	store at -20 °C
Note	For research use only
Application notes	SAP-Exo Kit is a quick, easy and reliable enzymatic cleanup reagent for PCR product cleanup. Removes excess primers and dNTPs; Fast 20 minute protocol; Add directly to PCR product; 100% Sample Recovery; Scalable for different reaction sizes; No interference on downstream applications; Easy to automate.
Expiration Date	12 months from date of receipt.