

Product Datasheet

Rab3AGST-His protein (orb533513)

Description	N-terminal tagged Rab3A is a small GTPase that belongs to the Ras superfamily. Rab proteins play an ...
Form/Appearance	liquid (Supplied in 50 mM Tris-HCl pH 8.0, 100 mM NaCl, 10 mM MgCl 2and 2 mM beta-mercaptoethanol)
Storage	store at -80 °C
Note	For research use only
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	E. coli
MW	Theoretical MW: 54 kDa
Expiration Date	6 months from date of receipt.