

Product Datasheet

Human Papillomavirus 6 GST protein (orb533531)

| | |
|------------------------|---|
| Description | Recombinant HPV6 antigen is a 55.6 kDa protein covering the full length of HPV6 major capsid, its N-... |
| Form/Appearance | liquid (Supplied in 50 mM arginine and 3M Urea) |
| Storage | store at -20 °C |
| Note | For research use only |
| Clonality | Recombinant |
| Purity | > 95% (SDS-PAGE) |
| Expiration Date | 6 months from date of receipt. |