

## Product Datasheet

### Human Papillomavirus 16 GST protein (orb533533)

<b>Description</b>	The Recombinant HPV-16 antigen is a full length protein expressed in E. coli having an MW of 56.7 kDa....
<b>Form/Appearance</b>	liquid (Supplied in PBS, 100 mM Arginine and 3 M Urea)
<b>Storage</b>	store at -20 °C
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Purity</b>	> 90% (SDS-PAGE)
<b>Expiration Date</b>	6 months from date of receipt.