

Product Datasheet

Rab9GST-His protein (orb533550)

| | |
|------------------------|---|
| Description | Rab9 is a small GTPase that belongs to the Ras superfamily. Rab proteins play an important role in v... |
| Form/Appearance | liquid (Supplied in 50 mM Tris-HCl pH 8.0, 100 mM NaCl, 10 mM MgCl ₂ , 10 μM GDP and 2 mM mercaptoethanol) |
| Storage | store at -80 °C |
| Note | For research use only |
| Clonality | Recombinant |
| Purity | > 90% (SDS-PAGE) |
| Source | E. coli |
| Expiration Date | 6 months from date of receipt. |