

Product Datasheet

PTP 1B protein (orb533679)

Description	PTP 1B is the prototype nontransmembrane protein tyrosine phosphatase. It is a ubiquitous protein t...
Form/Appearance	liquid (Supplied in 50 mM HEPES pH 7.4, 40 mM NaCl, 1 mM EDTA, 1 mM DTT and 2 mM beta-Mercaptoethanol)
Storage	store at -80 °C
Note	For research use only
Application notes	Activity: 1 nmol of PTP 1B will hydrolyze 37 nmol of pNPP per minute at 30°C.
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	E. coli
MW	Theoretical MW: 37.3 kDa
Expiration Date	6 months from date of receipt.