

Product Datasheet

FTaseW102T_Y154TGST protein (orb533684)

Description	FTase is one of the three mammalian protein prenyltransferases. The double mutant of FTase was engin...
Form/Appearance	liquid (Supplied in 25 mM HEPES pH 7.2, 40 mM NaCl and 2 mM DTT)
Storage	store at -80 °C
Note	For research use only
Application notes	Activity: 1 pmol of FTase double mutant will transfer 0.4 pmol of BGPP to H-Ras in 15 min at 37°C.
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	E. coli
MW	Theoretical MW: GST- α : 73 kDa, β : 48 kDa
Expiration Date	6 months from date of receipt.