

---

## Product Datasheet

### Parathyroid Hormone Antibody / N Terminal (orb534065)

**Description**

Parathyroid hormone (PTH), parathormone or



parathy...

**Species/Host**

Rabbit

**Reactivity**


Human

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P


 IHC testing  
of FFPE  
human  
parathyroid  
gl...

**Immunogen**

A synthetic peptide from the N-terminal region of human Parathyroid hormone was used as the immunogen for the recombinant Parathyroid Hormone antibody.

**Preservatives**

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

**Storage**

Store the recombinant Parathyroid Hormone antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

**Note**

For research use only

**Application notes**

The stated application concentrations are suggested starting points. Titration of the recombinant Parathyroid Hormone antibody may be required due to differences in protocols and secondary/substrate sensitivity.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

**Formula**

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

**Isotype**

Rabbit IgG, kappa

**Clonality**

Recombinant

**Purity**

Protein A affinity chromatography

**Uniprot ID**
**P01270**
**Hazard Information**

This recombinant Parathyroid Hormone antibody is available for research use only.

**Dilution Range**

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT, Prediluted IHC only format: incubate for 30 min at RT (1)

**Expiration Date**

12 months from date of receipt.