

Product Datasheet

OCT-2 Antibody / POU2F2 (orb534198)

Description

Transcription factor that specifically binds to the octamer


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

IHC testing of
FFPE human
tonsil with
OCT-2...

Immunogen

A portion of amino acids 112-297 from the human protein was used as the immunogen for the OCT-2 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the OCT-2 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).


Note

For research use only

Application notes

The stated application concentrations are suggested starting points. Titration of the OCT-2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Western blot
testing of
human
Ramos lysa...


Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity

Uniprot ID

P09086

Hazard
Information

This OCT-2 antibody is available for research use only.

Dilution Range

Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.