
Product Datasheet

RPA32 Antibody / Replication Protein A2 (orb534400)

Description

As part of the heterotrimeric replication protein A


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

WB

Western blot testing of human A375 cell ...

Immunogen

Full length recombinant human protein was used as the immunogen for this RPA32 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the RPA32 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The stated application concentrations are suggested starting points. Titration of the RPA32 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P15927

Hazard Information

This RPA32 antibody is available for research use only.

Dilution Range

Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.