
Product Datasheet

SLC2A1 antibody (orb534876)

Description

Rabbit monoclonal antibody to SLC2A1

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, FC, IHC-P

Immunogen

Recombinant fragment of human GLUT1 protein (around aa 203-305) (exact sequence is proprietary)

Target

SLC2A1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

 Antibody with azide - store at 2 to 8°C.
 Antibody without azide - store at -20 to -80°C.
 Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Rabbit / IgG

Clonality

Monoclonal

MW

55kDa

Uniprot ID
[P11166](#)
Entrez
[6513](#)
Dilution Range

ELISA: 1-2 µg/ml, IHC-P: 0.5-1 µg/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.