

## Product Datasheet

### MDM2 antibody (orb535761)

<b>Description</b>	Rabbit polyclonal antibody against MDM2 conjugated to HRP
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Synthetic peptide from the mid-protein of Human MDM2 (aa. 100-200)
<b>Target</b>	MDM2
<b>Preservatives</b>	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Conjugated antibodies should be stored at 4°C
<b>Note</b>	For research use only
<b>Application notes</b>	#REF!
<b>Clonality</b>	Polyclonal
<b>MW</b>	50 kDa
<b>Uniprot ID</b>	<a href="#">Q00987</a>
<b>NCBI</b>	<a href="#">NP_001138811.1</a>
<b>Entrez</b>	<a href="#">4193</a>
<b>Dilution Range</b>	WB (1:1000); IHC (1:50)
<b>Expiration Date</b>	12 months from date of receipt.