

Product Datasheet

BAS1 antibody (FITC) (orb53597)

Description	BAS1 antibody (FITC)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	FITC
Immunogen	Recombinant Arabidopsis thaliana 2-Cys peroxiredoxin BAS1, chloroplastic protein (66-266AA)
Target	BAS1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96291
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.