

## Product Datasheet

### Tau antibody (orb536193)

<b>Description</b>	Rabbit polyclonal antibody against Tau conjugated to HRP
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	Recombinant Human Tau 2NR4 P301S Fibril
<b>Target</b>	Tau
<b>Preservatives</b>	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Conjugated antibodies should be stored at 4°C
<b>Note</b>	For research use only
<b>Application notes</b>	A 1:1000 dilution of SPC-802 was sufficient for detection of Tau in 10 µg of human SH-SY5Y cell lysate cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
<b>Clonality</b>	Polyclonal
<b>MW</b>	45.8kDa
<b>Uniprot ID</b>	<a href="#">P10636</a>
<b>Entrez</b>	<a href="#">4137</a>
<b>Dilution Range</b>	WB (1:1000)
<b>Expiration Date</b>	12 months from date of receipt.