



## Product Datasheet

### NFKBIA antibody (orb536318)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to NFKBIA  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ICC, IF, WB   |
| <b>Immunogen</b>           | A synthetic peptide of human IκBα(NP_065390.1).   |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal,50% glycerol,pH7.3.   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 36kDa   |
| <b>Uniprot ID</b>          | <b>P25963</b>   |
| <b>Dilution Range</b>      | WB: 1:500-1:2000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Immunofluorescence analysis of NIH/3T3 c...



Immunofluorescence analysis of C6 cells ...