



Product Datasheet

ABL1 (Phospho-Y204) antibody (orb536407)


Description	Rabbit polyclonal antibody to ABL1 (Phospho-Y204)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y204 of human ABL1(NP_009297.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	122kDa/124kDa
Uniprot ID	P00519
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunofluorescence analysis of U2OS cell...



Immunofluorescence analysis of MCF-7 cel...