



Product Datasheet

DiMethyl-DNMT3A-K44 antibody (orb536463)


Description	Rabbit polyclonal antibody to DiMethyl-DNMT3A-K44
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	A synthetic dimethylated peptide around K44 of mouse DNMT3A(NP_031898.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	77kDa/101kDa
Uniprot ID	O88508
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of C6 cells ...