


Product Datasheet

Btk (Phospho-Y223) antibody (orb536505)

Description	Rabbit polyclonal antibody to Btk (Phospho-Y223)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y223 of human Btk (NP_000052.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q06187
Entrez	695
Dilution Range	WB: 1:500-1:2000, IF: 1:20-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...