



Product Datasheet

KRT7 antibody (orb536559)


Description	Rabbit polyclonal antibody to KRT7
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Cytokeratin 7 (KRT7)(NP_005547.3).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51KDa
Uniprot ID	P08729
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat lung...



Immunohistochemical staining of mouse lu...



Immunohistochemical staining of mouse ki...