



Product Datasheet

PRMT7 antibody (orb536560)


Description	Rabbit polyclonal antibody to PRMT7
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human PRMT7(NP_061896.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	62kDa/63kDa/73kDa/78kDa
Uniprot ID	Q9NVM4
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gallbladder tissue



Immunohistochemical staining of human stomach tissue



Immunohistochemical staining of human brain tissue