


Product Datasheet

BCL2 (Phospho-T56) antibody (orb536561)

Description	Rabbit polyclonal antibody to BCL2 (Phospho-T56)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around T56 of human Bcl-2(NP_000624.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22kDa/26kDa
Uniprot ID	P10415
Dilution Range	IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human br...