



## Product Datasheet

### ABL1 (Phospho-Y412) antibody (orb536562)


<b>Description</b>	Rabbit polyclonal antibody to ABL1 (Phospho-Y412)
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IF, IHC
<b>Immunogen</b>	A synthetic phosphorylated peptide around Y412 of human ABL1(NP_009297.2).
<b>Preservatives</b>	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	122kDa/124kDa
<b>Uniprot ID</b>	<b>P00519</b>
<b>Dilution Range</b>	IHC-P: 1:50-1:100, IF: 1:100-1:200
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence analysis of methanol-  
...



Immunofluorescence analysis of MCF-7  
cel...



Immunohistochemical staining of human  
br...