



Product Datasheet

RASSF1 antibody (orb536584)


Description	Rabbit polyclonal antibody to RASSF1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human RASSF1(NP_009113.3).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	8kDa/9kDa/16kDa/21kDa/31kDa/38kDa/39kDa
Uniprot ID	Q9NS23
Dilution Range	IHC-P: 1:50-1:100, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cel...



Immunohistochemical staining of mouse br...



Immunohistochemical staining of human li...