



Product Datasheet

ARFIP2 antibody (orb536601)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to ARFIP2 |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC |
| Immunogen | Recombinant fusion protein of human ARFIP2 (NP_036534.1). |
| Preservatives | Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Affinity purification |
| Uniprot ID | P53365 |
| Entrez | 23647 |
| Dilution Range | IHC-P: 1:50-1:200 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical staining of human brain tissue showing ARFIP2 expression.



Immunohistochemical staining of human lung tissue showing ARFIP2 expression.