

Product Datasheet

MTA2 antibody (orb536639)

Description Rabbit polyclonal antibody to MTA2

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen Recombinant fusion protein of human MTA2(NP_004730.2).

Preservatives PBS with 0.01% thiomersal,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 55kDa/75kDa


Uniprot ID **094776**

Dilution Range WB: 1:200-1:1000, IHC-P: 1:50-1:100


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of human tissue.