



Product Datasheet

ZNF433 antibody (orb536701)


Description	Rabbit polyclonal antibody to ZNF433
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	Recombinant fusion protein of human ZNF433(NP_001073880.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	73kDa/77kDa
Uniprot ID	Q8N7K0
Dilution Range	IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ut...



Immunohistochemical staining of human ga...



Immunohistochemical staining of human ap...