



Product Datasheet

ALOX12 antibody (orb536703)


Description	Rabbit polyclonal antibody to ALOX12
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human ALOX12(NP_000688.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	76KDa
Uniprot ID	P18054
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of rat lung ...



Immunohistochemical staining of mouse br...



Immunohistochemical staining of human co...