

Product Datasheet

ARHGEF5 antibody (orb536704)

Description Rabbit polyclonal antibody to ARHGEF5

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC

Immunogen Recombinant fusion protein of human ARHGEF5(NP_005426.2).

Preservatives PBS with 0.01% thiomersal,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 60kDa/176kDa


Uniprot ID [Q12774](#)

Dilution Range IHC-P: 1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of mouse kidney tissue.



Immunohistochemical staining of human tissue.



Immunohistochemical staining of human brain tissue.