



Product Datasheet

TYROBP antibody (orb536709)

Description	Rabbit polyclonal antibody to TYROBP
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human TYROBP(NP_003323.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11kDa/12kDa
Uniprot ID	O43914
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lung tissue.



Immunohistochemical staining of rat heart tissue.