



Product Datasheet

IL12RB1 antibody (orb536712)

Description	Rabbit polyclonal antibody to IL12RB1
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human IL12RB1(NP_714912.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42kDa/72kDa/73kDa
Uniprot ID	P42701
Dilution Range	IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF7 cell...



Immunohistochemical staining of rat kidn...