



Product Datasheet

UGDH antibody (orb536728)


Description	Rabbit polyclonal antibody to UGDH
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human UGDH(NP_003350.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44kDa/47kDa/55kDa
Uniprot ID	O60701
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse kidney



Immunohistochemical staining of human esophagus



Immunohistochemical staining of human lung