

Product Datasheet

IL36G antibody (orb536732)

Description Rabbit polyclonal antibody to IL36G

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Recombinant fusion protein of human IL36G(NP_062564.1).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 14kDa/18kDa


Uniprot ID [Q9NZH8](#)

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:100


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of mouse lu...



Immunohistochemical staining of human es...



Immunohistochemical staining of rat lung...