



Product Datasheet

MYC (Phospho-S62) antibody (orb536736)


Description	Rabbit polyclonal antibody to MYC (Phospho-S62)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S62 of humanMYC
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	62KDa
Uniprot ID	P01106
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lu...



Immunohistochemical staining of human ap...



Immunohistochemical staining of human lu...