

Product Datasheet

C12orf10 antibody (orb536752)

Description Rabbit polyclonal antibody to C12orf10

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC

Immunogen Recombinant fusion protein of human C12orf10(NP_067653.3).

Preservatives PBS with 0.01% thiomersal,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 42kDa


Uniprot ID **Q9HB07**

Dilution Range IHC-P: 1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human st...



Immunohistochemical staining of human br...



Immunohistochemical staining of human lu...