


Product Datasheet

PLCG1 (Phospho-Y783) antibody (orb536783)

Description	Rabbit polyclonal antibody to PLCG1 (Phospho-Y783)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y783 of human PLC gamma 1 (PLCG1)(NP_877963.1?).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	155KDa
Uniprot ID	P19174
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...