


Product Datasheet

JAK2 (Phospho-Y1007) antibody (orb536784)

Description	Rabbit polyclonal antibody to JAK2 (Phospho-Y1007)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Y1007 of human Phospho-JAK2-Y1007
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	120kDa
Uniprot ID	O60674
Entrez	3717
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...