

Product Datasheet

Bad (Phospho-S112) antibody (orb536803)

Description	Rabbit polyclonal antibody to Bad (Phospho-S112)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S112 of mouse Bad(NP_031548.1?).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18kDa
Uniprot ID	Q92934
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.