


Product Datasheet

CCND1 (Phospho-T286) antibody (orb536874)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to CCND1 (Phospho-T286) |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | A synthetic phosphorylated peptide around T286 of human Cyclin D1(NP_444284.1). |
| Preservatives | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 36KDa |
| Uniprot ID | P24385 |
| Dilution Range | WB: 1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of 293...