



Product Datasheet

CTNNB1 (Phospho-S33/S37/T41) antibody (orb536885)


Description	Rabbit polyclonal antibody to CTNNB1 (Phospho-S33/S37/T41)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A synthetic phosphorylated peptide around S33 & T41 of human beta Catenin(NP_001091679).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85kDa
Uniprot ID	P35222
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat skin...



Immunohistochemical staining of human lu...



Immunohistochemical staining of human co...