


Product Datasheet

P90RSK (Phospho-T359/S363) antibody (orb536940)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Description | Rabbit polyclonal antibody to P90RSK (Phospho-T359/S363) |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | A synthetic phosphorylated peptide around T359 & S363 of human p90Rsk/RSK1/RPS6KA1(NP_002944.2). |
| Preservatives | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 83kDa |
| Uniprot ID | Q15418 |
| Dilution Range | WB: 1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of A43...