


Product Datasheet

MET (Phospho-Y1234/1235) antibody (orb537031)

Description	Rabbit polyclonal antibody to MET (Phospho-Y1234/1235)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y1234 & Y1235 of human MET(NP_000236.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	145KDa
Uniprot ID	P08581
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HaL...