

Product Datasheet

SMAD3 (Phospho-T179) antibody (orb537041)

Description Rabbit polyclonal antibody to SMAD3 (Phospho-T179)

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen A synthetic phosphorylated peptide around T179 of human Smad3 (NP_005893.1).

Preservatives Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG

Clonality Polyclonal

Purity Affinity purification


MW 48kDa

Uniprot ID [P84022](#)


Entrez [4088](#)

Dilution Range WB: 1:500-1:2000


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of rat live...



Immunohistochemical staining of mouse ki...



Immunohistochemical staining of human co...