

Product Datasheet

Smad2 (Phospho-S465/467) + Smad3 (Phospho-S423/425) antibody (orb537116)

Description	Rabbit polyclonal antibody to Smad2 (Phospho-S465/467) + Smad3 (Phospho-S423/425)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic phosphorylated peptide around S465 & S467 of human Phospho-Smad2-S465/467 + Smad3-S423/425SMAD2(NP_005892.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60KDa
Uniprot ID	Q15796 , P84022
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.