


Product Datasheet

STAT1 (Phospho-Y701) antibody (orb537118)

Description	Rabbit polyclonal antibody to STAT1 (Phospho-Y701)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around Y701 of human STAT1(NP_009330.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	84KDa/91KDa
Uniprot ID	P42224
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa...