


Product Datasheet

Tyk2 (Phospho-Y1054/1055) antibody (orb537121)

Description	Rabbit polyclonal antibody to Tyk2 (Phospho-Y1054/1055)
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y1054 & Y1055 of human Tyk2(NP_003322.3).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	134KDa
Uniprot ID	P29597
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeL...