


Product Datasheet

TRIT1 antibody (orb537147)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to TRIT1 |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Recombinant fusion protein of human TRIT1(NP_060116.2). |
| Preservatives | PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 53kDa |
| Uniprot ID | Q9H3H1 |
| Dilution Range | WB: 1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of HeL...