



Product Datasheet

DSP antibody (orb537171)


Description	Rabbit polyclonal antibody to DSP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human Desmoplakin(NP_004406.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	315kDa
Uniprot ID	P15924
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse heart tissue



Immunohistochemical staining of human gastric tissue



Western blot analysis of extracts of HepG2 cells